



Applied Biological Materials Inc

Tel: 1-866-757-2414  
Email: info@abmGood.com  
Website: www.abmGood.com

## RNAitectin™

|          |             | Store at 4°C |
|----------|-------------|--------------|
| Cat. No. | Description | Quantity     |
| G073     | RNAitectin  | 1.0ml        |

### Description

RNAitectin™ is a transfection reagent specially formulated with multiple cationic polymers. It is suitable for the transfection of RNAi oligos into cultured eukaryotic cells.

### Transfection protocol

Use the following conditions as guidelines to transfect mammalian cells in a 6-well or 35mm dish format. For other culture vessels, please refer to Table 1.

1. Eighteen to twenty-four hours prior to transfection, seed cells at a density of  $1-3 \times 10^5$  cells per well in 2.0ml of appropriate growth medium (with serum and antibiotics if normally required). Incubate the cells at 37°C in a CO<sub>2</sub> incubator until cells are 70% to 90% confluent at the time of transfection.

**Suspension Cells:** Just prior to preparing complexes, plate  $3-5 \times 10^5$  cells in 0.8ml of serum-free medium without antibiotics.

Since transfection efficiency is sensitive to culture confluence, it is important to maintain a standard seeding protocol throughout all of your experiments.

2. For each transfection sample, prepare the complexes in sterile micro centrifuge tubes as follows:

**Solution A:** Dilute 1-3µg of RNAi oligos into 125µl of serum-free, antibiotic-free medium.

**Solution B:** Mix RNAitectin™ reagent thoroughly prior to use, then dilute 4-10µl of reagent in 125µl serum-free, antibiotic-free medium.

Incubate solution A and B separately at room temperature for 5 minutes.

3. Combine solutions A and B, mix thoroughly to ensure uniform distribution and incubate for 20 minutes at room temperature to allow RNA-liposome complexes to form.

4. **Adherent cells ONLY (For suspension cells, go directly to step 5b):**

Add 0.8ml of serum-free, antibiotic-free medium to the RNAitectin™-RNA complex. Mix solution gently.

- 5a. **Adherent cells:** Remove growth medium from the cells and add 1.0ml of RNAitectin™-RNA solution to the each well containing cells. Incubate at 37°C.

- 5b. **Suspension cells:** Add 0.2ml of the RNAiVectin™ -RNA solution into each well containing suspension cells in 0.8ml serum-free, antibiotic-free medium. Incubate at 37°C.
- After 5-8 hours, remove transfection solution and add 2.0ml of the appropriate growth medium (with serum and antibiotics) or add 0.1ml of FBS directly into each vessel. Incubate the cells at 37°C for a total of 18-24 hours.
  - Assay cell extracts for marker gene activity 24-72 hours after the start of transfection depending on cell type and promoter activity.
  - A similar procedure can be used to transfect an RNAi vector for stable expression. Seventy-two hours after transfection passage the cells 1:10 or higher into an appropriate selective medium.

### Optimizing transfection for specific cell lines

To achieve the maximum transfection efficiency and low cytotoxicity, optimize transfection conditions by varying cell density along with RNA oligo and RNAiVectin™ concentrations. Optimal results have been observed when cells are 70-90% confluent and RNA(µg): RNAiVectin™ (µl) ratios are 1:1 to 1:5.

**Table 1: Reagent quantities for different culture vessels**

| Culture Vessel | Surface area per well (cm <sup>2</sup> ) | Volume of plating medium | RNAi oligos (µg) in medium volume (µl) | DNAiVectin™ in medium volume | Transfection medium vol. |
|----------------|------------------------------------------|--------------------------|----------------------------------------|------------------------------|--------------------------|
| 24-well        | 2                                        | 500 µl                   | 0.2-0.4 µg in 25 µl                    | 2-4 µl in 25 µl              | 0.4 ml                   |
| 12-well        | 4                                        | 1 ml                     | 0.5-0.8 µg in 100 µl                   | 5-8 µl in 100 µl             | 0.6 ml                   |
| 6-well         | 10                                       | 2 ml                     | 1.0-3.0 µg in 125 µl                   | 4-10 µl in 125 µl            | 0.8 ml                   |
| 35mm           | 10                                       | 2 ml                     | 1.0-3.0 µg in 125 µl                   | 4-10 µl in 125 µl            | 0.8 ml                   |
| 60mm           | 20                                       | 5 ml                     | 3.0-6.0 µg in 500 µl                   | 12-30 µl in 500 µl           | 2.4 ml                   |
| 100mm          | 60                                       | 10 ml                    | 8.0-16.0 µg in 1.5 ml                  | 32-80 µl in 1.5 ml           | 6.4 ml                   |

### Notes

- Do not add antibacterial agents to the media during transfection
- Cells that will not tolerate the absence of serum for 2-24 hours can be transfected in the presence of serum. This is done by preparing the oligo-liposome complexes for 45 minutes in serum-free medium, followed by diluting the complexes with serum-containing medium before adding to the cells. It is extremely important that the amount of RNAiVectin™ be re-optimized as the optimal amount of lipid under these conditions may be different from that observed for serum-free transfection.

*For laboratory research only. Not for clinical applications.  
For technical questions, phone the ABM helpline at 1-866-757-2414  
or visit our website at [www.abmGood.com](http://www.abmGood.com)*

**CERTIFICATE OF ANALYSIS**