



image not found or type unknown

**Applied Biological Materials Inc.**

Telephone: 1-866-757-2414

Email: [info@abmgood.com](mailto:info@abmgood.com)

Website: [www.abmGood.com](http://www.abmGood.com)

---

## Universal 3' miRNA Reverse Primer

**Cat. No.**  
MPH00000

**Unit**  
150 µl / 10 µM

---

**Cat. No.** MPH00000

**Name** Universal 3' miRNA Reverse Primer

**Unit** 150 µl / 10 µM

**Description** The Universal 3' miRNA Reverse Primer is compatible with miRNA cDNA synthesized using **abm's** [miRNA All-In-One cDNA Synthesis Kit](#) (Cat. No. G898). This reverse primer is specific for the unique and proprietary sequence incorporated into the miRNA cDNA by the Oligo d(T) with Adapter in the [miRNA All-In-One cDNA Synthesis Kit](#). miRNA cDNA lacking this adapter sequence will not be compatible with this reverse primer.

This product is provided at 10 µM and can be paired with any of **abm's** [miRNA qPCR forward primers](#) for qPCR expression analysis using [Blastaq™ 2X qPCR Mastermix](#) (Cat. No. G891). Please note that this reverse primer is also provided complimentary with every miRNA forward primer ordered.

A variety of housekeeping gene primers for reaction controls are also available including: U6-2 (Cat. No. [MPH00001](#)), SNORD44 (Cat. No. [MPH00003](#)), and SNORD48 (Cat. No. [MPH00005](#)).

**System** qPCR Primers

**Storage Condition** -20°C

---

**Caution:** *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

*Please contact a technical service representative for more information (1-866-757-2414).*