



image not found or type unknown

Applied Biological Materials Inc.

Telephone: 1-866-757-2414

Email: info@abmgood.com

Website: www.abmGood.com

U6 Primers

Cat.No.	Unit	Accession Number
MPH00001	150 µl / 10 µM	NR_002752, NR_003027, X07425.1

Cat. No.	MPH00001
Name	U6 Primers
Unit	150 µl / 10 µM
Description	abm miRNA primers are for qPCR applications and each miRNA-specific primer set comes with a universal reverse primer. All miRNA primers are supplied in 10µM format.
miRNA Name	U6-2, U6
Accession Number	NR_002752, NR_003027, X07425.1
Species	Human, Mouse
Note	abm's miRNA All-In-One cDNA Synthesis Kit (Cat No. G898) must be used to reverse transcribe the RNA. Failure to do so will not yield valid qPCR results. Each primer set includes a universal reverse primer and is specific for the unique and proprietary sequence incorporated into the cDNA by the miRNA adapter in the miRNA cDNA synthesis kit. The miRNA qPCR primer sets cannot accurately detect cDNA generated using other first strand synthesis kits.

Caution: *This product is for research use only and is not intended for therapeutic or diagnostic applications.*

Please contact a technical service representative for more information (1-866-757-2414).